



HM-B19, 3330 Hospital Drive NW
Calgary, Alberta Canada T2N 4N1
T: (403) 220.5511
F: (403) 270.2964
E: acft.lab@ucalgary.ca

APPENDIX VI

URINE RECOLLECTION

This document is for collection sites to use if urine recollection is required if the following situations occurred at initial attempt of collection.

- Urine temperature is outside the acceptable range.
- Urine volume is lower than 30 mL.
- Client is unable to void.

Ideally, urine recollection should occur on the same day or next day. Provide another labelled sample container to client. Give the client a glass of water, wait 30 minutes, and have the client try again.

If the client cannot void during the initial attempt and refuses to provide a recollection, send the empty specimen cup to ACFT. Make a Lab Comment in Epic or on the ACFT requisition.

If the client is unable to void at the second attempt (recollection) on the same day or the following day, send the initial sample cup to ACFT. Make a Lab Comment in Epic or on the ACFT requisition.

If samples were recollected due to temperature or volume issues, a single registration (or one ACFT requisition) is adequate for the recollection.

Indicate on the specimen container which is the initial collection, and which is the recollection prior to sending the samples to ACFT.

For questions regarding the urine collection please contact ACFT. For questions regarding Epic please contact the AHS Epic Team.