

DRY POWDER INHALER

ELLIPTA



1 SLIDE OPEN AND HEAR "CLICK"



2 BREATHE OUT




3 DEEP BREATH IN AND HOLD



4 DO NOT BLOCK AIR VENTS



5 BREATHE OUT



6 CLOSE INHALER



7 RINSE, GARGLE AND SPIT

IMPORTANT INFORMATION

- Note:** Place mouthpiece between your teeth and close your lips.
- Empty?:** The window has numbers that count down the dose – "0" means it is empty. Throw it away.
- Expired?:** Check the expiry date on the back of the device.
- Reminders:** Keep the device dry.
Do not breathe into the device.
Take once daily – same time every day.
Do not block air vents with your fingers.
- Cleaning:** Wipe mouthpiece with a dry tissue or cloth.

This material is designed for information purposes only. It should not be used in place of medical advice, instruction and/or treatment. If you have specific questions, please consult your doctor or appropriate healthcare professional.

Illustrations developed by the Community Pediatric Asthma Service, www.ucalgary.ca/icancontrolasthma

For more information on inhaled medications and device demos, visit our website at:
www.ucalgary.ca/asthma



Supporting current
Respiratory Guidelines

